


526 – Metal Detectable Broom/Squeegee Handle

Document Number: HSE 7237

Date Of Issue: 26/07/23

Revision Number: 1

Date of Revision: 26/07/23



Code	Description	Colour	Size (mm)	Weight (kg)
510-P01-A43	Ultra Hygienic Blue Plastic Broom/Squeegee Handle Only	Blue	1500 ø 32	0.520

Strong and resistant one piece handle, the knob and the 32 mm diameter handle make it ergonomic and easy to use. Equipped with 22 mm Italian thread. Tests show a deformation of the handle under pressure of less than 60%.

Minimum Temperature -20°C

Maximum Temperature 100°C

Autoclave 121°C for 20mins


User and Maintenance Manual

- Remove packaging and any label from tools and disenfect them
- Both manual and autoclave disinfection is possible
- Sanitize tools using appropriate chemical decontamination products, in the concentrations, times and temperatures provided
- Make sure tools completely dry up before next use
- Hang tools while drying up to avoid any contamination or breakage
- The storing area must be clean and far away from sources of intense light or heat
- Do not use damaged tools

Compliance

- REACH
- FDA
- EU 1935/2004 and updates
- EU 2023/2006 and updates
- EU 10/2011 and updates

No warranty is given or implied with respect to this information or patent infringement. Detectamet Ltd do not accept liability for loss or damage arising from the use of this information. Results are based on a test sample, our general experience and information from our suppliers. Data and results must be confirmed by the buyer by testing for its intended conditions of use.



Helen Morrison
Group Managing Director